

Make Alarm Configuration
easy. 

07

Final design

This chapter presents the final design: SetWise. Initially, an overview of the final design is provided, followed by a detailed walkthrough of each interface screen, including the key features. The chapter concludes with an explanation of how SetWise fits into the broader system architecture, including an explanation of the underlying algorithms.

7.1 The Design: SetWise

7.2 System architecture




SetWise


[Click here to interact](#)

SetWise is a **configuration interface** designed to support nurses in creating a personalised set of alarm rules for each patient, intending to reduce the overall alarm load in the ICU. While simplifying the configuration process, SetWise generates alarms that are both clinically relevant and context-aware. The interface acts as input for multiple algorithms to create the personalised alarm rules. SetWise is located next to the patient bed and can also be accessed from computers at the central nurse station.

Two-Step Configuration

SetWise simplifies the setup process by guiding nurses through just two essential steps. Select the patient's primary clinical problem and risk level, and the system automatically configures detailed monitoring settings.



The screenshot illustrates the SetWise configuration interface, which consists of two main sections:

- Top Section (Patient Overview):** This section includes the patient's name (S. Groen), gender (Vrouw), age (16 jaar), and weight (55 kg). It also shows the primary clinical problem (Respiratoire Insufficiëntie) and risk level (Mid). Below this, additional information is provided, such as Sepsis, Diabetes, Herstellende, Pneumonie, and Chronisch Nierfalen.
- Bottom Section (Monitoring op Orgaan Niveau):** This section displays monitoring levels for various organ systems: Circulatoir (mid), Respiratoir (strik), Neuro (los), Lever (los), Nieren (los), and Overig (mid). Below each organ system, specific monitoring parameters are listed:
 - Circulatoir: HR 70 - 110, BP Mean 65 - 85
 - Respiratoir: Saturatie 90 - 100, AF 10 - 25
 - Neuro: ICP NIRS
 - Lever: ALAT ASAT
 - Nieren: KREAT Ureum
 - Overig: Temp. 36.5 - 39
 Each parameter has an "Aanpassen" (Adjust) button below it.

Clear overview

SetWise provides a comprehensive overview of the configured monitoring settings across all organ systems. This layout enables nurses to quickly verify or adjust settings without losing sight of the bigger clinical picture.

Context-Aware Shortcuts

Use the predefined additional information labels to add contextual details about the patient easily. These inputs further refine the configuration of target ranges and detailed monitoring levels, ensuring that the alarm settings accurately reflect the patient's current situation.

Dive deep, if desired

The proposed detailed monitoring settings can always be fine-tuned.

Patient Overview - Home

The patient overview serves as the landing page of the interface. It is the starting point for configuring a new patient, and acts as the home page to return to for an overview of all patients in the unit.

Patiënten Overzicht

The home page enables an overview of all patients and their configured risk level, indicated by the colour of the beds.


Risk Level Legend:

- High Risk Level (Dark Blue)
- Mid Risk Level (Medium Blue)
- Low Risk Level (Light Blue)

Add a new patient
To add a new patient, click the plus button of an empty bed.

Patient Specific Settings
Easily navigate to the alarm settings of a specific patient by clicking the 'Alarm Instellingen' button.

Patient information
View detailed patient information by clicking on the patient bed. This contextual information helps nurses make informed decisions about configuration.



Automated Patient Proposal
Based on data from the Electronic Health Record, the system automatically proposes the patient whose admission is next. The nurse can confirm and start configuration, or choose another patient from the list.

Configure a New Patient

New patients are configured in just two mandatory steps. As the user progresses through the configuration, the detailed organ-system monitoring settings adjust in real-time, allowing for immediate insight into the **implications of each configuration choice**.

Step-by-step Guidance

The progress bar guides the user step-by-step through the configuration process.

The Higher the Risk, the More Alarms

The second step is to assign a risk level to the patient. This classification is used by the system's underlying algorithms to determine the sensitivity and frequency of alarm triggers.

The main problem

The first step is to select the patient's main clinical problem, which reflects the **main** reason for admission. This can always be updated in the general alarm settings overview if the patient's condition changes during the ICU stay.

Instellen Nieuwe Patiënt

Wat is het hoofdprobleem?
Respiratoire Insufficiëntie

Wat is het risico - level van de patiënt?
 Laag Mid Hoog
sneller een alarm

Selecteer eventueel additionele informatie:
 Sepsis
 Diabetes
 Herstellende
 Pneumonie
 Chronisch Nierfalen

Monitoring op Orgaan Niveau

Circulatoir	Respiratoir	Neuro	Lever	Nieren	Overig
 mid	 strik	 los	 los	 los	 mid
HR 70 - 110 BP Mean 65 - 85	Saturatie 90 - 100 AF 10 - 25	ICP NIRS	ALAT ASAT	KREAT Ureum	Temp. 36.5 - 39

Terug Bevestigen

Real-Time Adjustment

Once the main problem and risk level are selected, the system automatically adjusts the monitoring levels and corresponding target ranges for each organ system in real-time.

Target Ranges

Each parameter displays a target range, which represents the values considered "normal" for a specific patient's condition, based on the main problem, additional information, and patient details.

Additional information

ICU patients can have complex pathologies and medical backgrounds. Nurses can further refine the configuration by selecting additional clinical information labels, and the detailed monitoring will automatically adjust.

Additional Information Labels

After a new patient is configured, the user is redirected to the patient overview screen. From there, the alarm settings can be accessed and refined, such as the Additional Information Labels. These features enable users to add detailed clinical information, allowing for more precise configuration of individual patients.

Search all
Search across all categories using the general search bar.

Filter on organ system
The categories help to find the right information.

Selection on top
Selected labels are shown at the top of the screen, offering a clear overview of all applied labels.

Suggested Labels
Because the database of available labels is extensive, the interface shows only a limited set of suggestions based on the selected clinical problem and patient background. Users can expand this list by opening the pop-up.

Real-time Feedback
As additional labels are added, the system updates the monitoring levels and target ranges immediately to show the impact of each label.

Circulatory Settings

By going one level deeper into the settings, users can adjust the target ranges of vital parameters for the circulatory organ system.

Never out of sight

The patient's information remains visible at the top of the screen to ensure that changes are always made for the right patient.

Patiënten Overzicht > Alarm Instellingen - S. Groen > Circulatoire Monitoring

Circulatoire Monitoring

Fine Tune Organ Monitoring

Although the system initially determines the monitoring level based on the patient's overall risk and clinical context, nurses can still manually adjust this level (loose, mid, or tight). These levels determine the system's sensitivity to deviations in vital parameters, which is visualised by the real-time adjustment of the target range graphs.

Monitoring Level

Los Mid Strak

Circulatoire Streefgebieden Aanpassen

HR

130 110 70 40

0 5 10 15 20 25 30 Secondes

BP Mean

100 85 65

0 5 10 15 20 25 30 Secondes

Geen alarm Geel alarm Rood alarm

Circulatoire Additionele informatie

Sepsis Hartfalen Bradycardie ...

Alarm Instellingen ←



Naar Respiratoir ↗

Tailored Additional Information

Only relevant additional information labels that affect circulatory monitoring are shown.

Next Organ System

Directly navigate to the other organ-systems.



No alarm, within safe range

Medium priority alarm

Critical Alarm

Target Range Graphs

The graphs show how long it takes before a vital parameter exceedance will trigger an alarm. For example, if the monitoring level is set to loose, a small deviation will trigger an alarm after approximately 15 seconds. **Note:** if the deviation is large, the alarm may be triggered sooner.

[Patiënten Overzicht > Alarm Instellingen - S. Groen > Circulatoire Monitoring](#)

Circulatoire Monitoring

Monitoring Level: Strak

Circulatoire Additionele informatie: Sepsis, Hartfalen, Bradycardie, ...

Circulatoire Streefgebieden Aanpassen

HR: 130 (Geen alarm), 110 (Geel alarm), 70 (Rood alarm)

BP Mean: 100 (Geen alarm), 85 (Geel alarm), 65 (Rood alarm)

See what you Set
Selecting an additional information label instantly updates **both** the target ranges AND the monitoring level, to show the implications of the label.

Automatically adjusted

Naar Respiratoir

Respiratory & Other Settings

The Respiratory and Other settings follow the same structure as the circulatory settings to ensure a streamlined experience.

Keep track of where you are

The “breadcrumb” navigation helps to stay oriented and return to the patient or alarm overview with ease.

[Patienten Overzicht > Alarm Instellingen - S. Groen > Overige Monitoring](#)

Overige Monitoring

Monitoring Level: Los, Mid, Strak (Strak is selected)

Overige Streefgebieden Aanpassen

Temperatuur

Waarde	Geen alarm	Geel alarm	Rood alarm
39	Green	Yellow	Red
36.5	Green	Yellow	Red

0 5 10 15 20 25 30 Secundes

34567 - S. Groen - Vrouw - 16 jaar - 55 kg
Respiratoire Insufficiëntie - Risico: mid

Overige Additionele informatie

- Delier
- Diabetes
- Herstellende
- ...

[Patienten Overzicht > Alarm Instellingen - S. Groen > Respiratoire Monitoring](#)

Respiratoire Monitoring

Monitoring Level: Los, Mid, Strak (Strak is selected)

Respiratoire Streefgebieden Aanpassen

Saturatie

Waarde	Geen alarm	Geel alarm	Rood alarm
100	Green	Yellow	Red
90	Green	Yellow	Red

0 5 10 15 20 25 30 Secundes

AF

Waarde	Geen alarm	Geel alarm	Rood alarm
40	Green	Yellow	Red
25	Green	Yellow	Red
10	Green	Yellow	Red

0 5 10 15 20 25 30 Secundes

34567 - S. Groen - Vrouw - 16 jaar - 55 kg
Respiratoire Insufficiëntie - Risico: mid

Respiratoire Additionele informatie

- Pneumonie
- Astma
- COPD
- ...

Alarm Instellingen ←


Naar Overig ⏪

Hover & Slide

By hovering over the slider, arrows appear to guide the user to adjust to the desired target range.

7.2 System Architecture

The interface is part of a larger ecosystem, namely the envisioned future alarm system of UMC Utrecht (see Section 3.5.2). Figure 7.1 illustrates how the interface integrates with the other components. The system is not yet ready for deployment; the Smart Alarm Manager software still needs to be developed by Ascom, and the Brain module is still several years from completion. Nonetheless, this user interface enables UMC Utrecht to test the envisioned system with nurses, which is essential for deployment.



Component Functionality

The brain automatically populates the configuration interface with the patient's contextual information, the nurse only needs to verify and/or adjust this information. Using this context, the cerebellum then applies the appropriate alarm protocols; for example, it will allow a lower heart rate for patients with bradycardia. **The nurse stays in control, and may always fine tune** these target ranges. For an initial deployment, the system could also function without the Brain component. The configuration interface can directly request the Cerebellum to select the right alarm protocol, using the main problem and additional information provided by the nurse. In either scenario, the contextual data, configuration settings, and chosen protocol are forwarded to the Smart Alarm Manager software, where the core algorithms are executed.

The primary purpose of these algorithms is to generate a set of patient-specific alarm rules, which the Smart Alarm Manager then uses to determine whether an incoming alarm should actually be raised. The algorithms use the input of the configuration interface (and Brain) to create these alarm rules. The depth-duration algorithm combines contextual information with the risk level of organ systems to calculate both how deep and for how long each vital parameter may deviate from its target. Meanwhile, the data-knowledge platform utilises alarm protocols and contextual patient information (including the main problem and additional labels) to generate personalised target ranges. Together, these algorithms produce the combined output: personalised alarm rules to eventually reduce the alarm burden in the ICU.

Accessible anywhere, from any screen

[Go to Prototype](#)